



Slite Tech Instruments




GUMMIKOFFERKLAMMERZANGE
KATALOG




SLITE TECH INSTRUMENTS GERMANY

RUBBER DAM CLAMP FORCEPS




CLAMPS (PLIERS)





RUBBER DAM CLAMP FORCEPS



2510 - 9
Winged
Labial Clamp



2511 - 9T
Winged
Edged Labial Clamp



2512 - 212
Winged
Cervical Clamp



2513 - 210
Winged
Labial Clamp

Premolars



2514 - 00
Winged



2515 - 0
Winged



2516 - W0
Wingless

For Small Teeth, Elongated or With Narrow Collars



2517 - 1
Winged




2518 - 1T
Winged



2519 - W1
Wingless

For General Use With Adjustable Jaws According
To Gingival Line Shapes



2520 - 2A
Winged



2521 - 2AT
Edged, Winged



2522 - W2A
Wingless

For Premolars, With Elongated Collar



2523 - 2
Winged



2524 - 2T
Edged, Winged



2525 - W2
Wingless

For Premolars, With Elongated Collar



2526 - 29
Wingless




2527 - 60
Wingless



RUBBER DAM CLAMP FORCEPS



2528 - 3
Wingless



2529 - W3
Wingless



2530 - 4
Winged



2531 - W4
Wingless



For Small Lower Molars



2532 - 7
Winged



2534 - 8
Winged



2535 - W8
Wingless


For Small Lower Molars With Round clamps

For Square Lower Molars


For Round Upper Molars



2536 - 8A
Winged



2537 - W8A
Wingless



2538 - 12A
Winged



2539 - 13A
Wingless

For Partially Broken Molars

For Molar Roots



2540 - 14
Winged



2541 - 14T
Edged Winged



2542 - W14
Wingless



2543 - 14A
Winged



2544 - W14A
Wingless

For Round Molars, Partially Broken, Undersized or Inclined



2545 - 26
Wingless




2246 - 51
Wingless



2547 - 202
Winged



2548 - 205
Winged





KONTAKTIERE JETZT



info@slitetechnstruments.de



+49 162 9224643



slitetechnstruments01@gmail.com



Osterkamp 10 33605,
Bielefeld Germany



www.slitetechnstruments.de