



Slite Tech Instruments



**NAHTINSTRUMENTE**  
**KATALOG**



# SLITE TECH INSTRUMENTS GERMANY

## TABLE OF CONTENT

01	Needle Holders	03-05
----	----------------	-------

02	Needle Holders with Tc	06-12
----	------------------------	-------

03	Suture Instruments	13
----	--------------------	----

04	Ligature Needles	14-18
----	------------------	-------

05	Ligature Conductors	19-22
----	---------------------	-------



# SLITE TECH INSTRUMENTS GERMANY


## NEEDLE HOLDERS

AF



12,0cm  
AF.300.12

12,0cm  
AF.301.12



12,5cm  
AF.302.12




13,0cm  
AF.304.13

13,0cm  
AF.305.13



14,0cm  
AF.314.14





18,0cm  
AF.324.18



# SLITE TECH INSTRUMENTS GERMANY

## NEEDLE HOLDERS



AF





## NEEDLE HOLDERS

AF



### Mathieu

Article No.	Size
AF.380.14	14,0cm
AF.380.17	17,0cm
AF.380.20	20,0cm



### Mathieu grazil

Article No.	Size
AF.380.84	14,0cm



### Castroviejo gerade

Article No.	Size
AF.430.14	14,0cm
AF.430.18	18,0cm
AF.430.22	22,5cm



### Castroviejo gebogen

Article No.	Size
AF.431.14	14,0cm
AF.431.18	18,0cm
AF.431.22	22,5cm



# SLITE TECH INSTRUMENTS GERMANY

## NEEDLE HOLDERS WITH TC



Mit Hartmetall-einsätzen  
*With Tungsten Carbide inserts*  
Avec mors en métal dur  
Con puntas de metal duro  
Con punte in metallo duro



### Crile-Wood

Article No.	Size
AF.762.15	15,0cm
AF.762.18	18,0cm
AF.762.20	20,0cm
AF.762.23	23,0cm



AF



# SLITE TECH INSTRUMENTS GERMANY

## NEEDLE HOLDERS WITH TC

AF




**Baby-Crile-Wood**

Article No.	Size
AF.763.15	15,0cm




**Schwed. Modell**

Article No.	Size
AF.766.15	15,0cm



**Sarot**

Article No.	Size
AF.766.18	18,0cm
AF.766.20	20,0cm
AF.766.26	26,0cm



**Mayo-Hegar**

Article No.	Size
AF.768.14	14,0cm
AF.768.16	16,0cm
AF.768.18	18,0cm
AF.768.20	20,0cm
AF.768.24	24,0cm
AF.768.26	26,0cm



# SLITE TECH INSTRUMENTS GERMANY

## NEEDLE HOLDERS WITH TC



Mayo-Hegar 20cm  
starke Ausführung

Article No.	Size
AF.770.20	20,0cm



MEDE - Vascular

15,0cm  
AF.771.15

18,0cm  
AF.771.18

20,0cm  
AF.771.20

23,0cm  
AF.771.23



De Bakey

18,0cm  
AF.772.18

20,0cm  
AF.772.20

23,0cm  
AF.772.23

26,0cm  
AF.772.26

31,0cm  
AF.772.31

**AF**


Mit Hartmetall-  
einsätzen  
*With Tungsten  
Carbid inserts*  
Avec mors en métal dur  
Con puntas de metal duro  
Con punte in metallo duro



# SLITE TECH INSTRUMENTS GERMANY


## NEEDLE HOLDERS WITH TC

AF




**Masson**

Article No.	Size
AF.776.27	27,0cm



**Micro-Ryder**

Article No.	Size
AF.782.13	13,0cm
AF.782.15	15,0cm
AF.782.18	18,0cm
AF.782.20	20,0cm
AF.782.23	23,0cm




**Ryder (Frenche eye)**

Article No.	Size
AF.780.15	15,0cm
AF.780.18	18,0cm
AF.780.20	20,0cm

glattes Maul  
*smooth jaws*  
bocas lisas  
à mors lisses


**Ryder (Frenche eye)**

Article No.	Size
AF.781.18	18,0cm



**Wangensteen**

Article No.	Size
AF.790.27	27,0cm



**Finochietto**


Article No.	Size
AF.791.27	27,0cm



# SLITE TECH INSTRUMENTS GERMANY

## NEEDLE HOLDERS WITH TC

AF



### Bozemann

Article No.	Size
AF.796.20	20,0cm
AF.796.24	24,0cm
AF.796.26	26,0cm



Mit Hartmetall-  
einsätzen

With Tungsten  
Carbid inserts

Avec mors en métal dur  
Con puntas de metal duro  
Con punte in metallo duro



### Toennis TC

Article No.	Size
AF.820.18	18,0cm



### Langenbeck TC

Article No.	Size
AF.824.16	16,0cm
AF.824.18	18,0cm
AF.824.20	20,0cm
AF.824.23	23,0cm



# SLITE TECH INSTRUMENTS GERMANY

## NEEDLE HOLDERS WITH TC

AF



### Olsen-Hegar

Article No.	Size
AF.800.12	12,0cm
AF.800.14	14,0cm
AF.800.17	17,0cm
AF.800.19	19,0cm



Mit Hartmetall-  
einsätzen

*With Tungsten  
Carbid inserts*

Avec mors en métal dur  
Con puntas de metal duro  
Con punte in metallo duro



### Mini-Gillis

Article No.	Size
AF.802.12	12,0cm

### Gillis

Article No.	Size
AF.802.15	15,0cm





# SLITE TECH INSTRUMENTS GERMANY

## NEEDLE HOLDERS WITH TC

AF



### Mathieu

Article No.	Size
AF.812.14	14,0cm
AF.812.17	17,0cm
AF.812.20	20,0cm



### Mathieu

Article No.	Size
AF.814.17	17,0cm
AF.814.20	20,0cm






# SLITE TECH INSTRUMENTS GERMANY

## SUTURE INSTRUMENTS

AF



**Hegenbarth**

Article No.	Size
AF.886.12	12,5cm



**Collin**

Article No.	Size
AF.892.13	13,0cm



**Michel**

Article No.	Size
AF.890.12	12,0cm



**Michel**

Article No.	Size
AF.891.13	13,0cm




**Klammerentferner**

Article No.	Size
AF.896.12	12,0cm

### FUNKTIONSWEISE | HOW TO USE

Schieben Sie vorsichtig den untern Führungsarm des Klammerentferners unter die Klammer.



Insert the nose of the Manipler skin staple Remover firmly under the staple.

Schließen Sie den K-Entferner so das die Klammer anfängt sich zu verformen.



To reform the staple for removal.

Der K-Entferner muß vollkommen geschlossen sein so das die Klammer ganz geöffnet ist.



completely squeeze the handles until they are fully closed.

Erst jetzt können Sie die den Entferner samt Klammer von der Stelle entfernen.



With the handles closed, lift the staple straight up from the tissue.



# SLITE TECH INSTRUMENTS GERMANY


## LIGATURE NEEDLES

AF



Schmieden-Dick

Article No.	Size
AF.920.32	32,0cm



Schmieden

Article No.	Size
AF.922.21	21,0cm

Biegsam  
*malleable*  
maleable  
malléable  
malleabile



Article No.	Size
AF.924.20	20,0cm



Copper

Article No.	Size
AF.926.19	19,0cm



## LIGATURE NEEDLES

AF



Stumpf  
*Blunt*  
Smusse  
Romas  
Mousses




Kronecker

Article No.	Size
AF.940.13	13,0cm

Kronecker

Article No.	Size
AF.942.13	13,0cm



Stumpf  
*Blunt*  
Smusse  
Romas  
Mousses

Deschamps

Article No.	Size
AF.950.24	24,0cm

Deschamps

Article No.	Size
AF.952.24	24,0cm



# SLITE TECH INSTRUMENTS GERMANY

## LIGATURE NEEDLES

AF



**Stumpf**  
**Blunt**  
Smusse  
Romas  
Mousses



**Deschamps**  
Article No.  
AF.956.20 | Size  
20,0cm



**Scharf**  
**sharp**  
agudo  
pointue  
acuti



**Deschamps**  
Article No.  
AF.957.20 | Size  
20,0cm


**Deschamps**  
Article No.  
AF.958.20 | Size  
20,0cm



# SLITE TECH INSTRUMENTS GERMANY

## LIGATURE NEEDLES

AF



Deschamps

Article No.

AF.967.28

Size

28,0cm

Stumpf

Blunt

Smusse

Romas

Mousses

Deschamps

Article No.




AF.968.28

Size

28,0cm



## LIGATURE NEEDLES




AF



# SLITE TECH INSTRUMENTS GERMANY

## LIGATURE CONDUCTORS

AF



Payer



Payer



# SLITE TECH INSTRUMENTS GERMANY

## LIGATURE CONDUCTORS



**Brunner Hohlsonde**

Article No.	Size
AF.987.18	18,0cm



**Kirschner Hohlsonde**

Article No.	Size
AF.988.27	27,0cm



9mm breit

**Brunner Hohlsonde**

Article No.	Size
AF.987.27	27,0cm

AF





# SLITE TECH INSTRUMENTS GERMANY

## LIGATURE CONDUCTORS

AF



**Kocher**

Article No.  
**AF.986.01**



Size  
**16,0cm**



**Kocher**

Article No.  
**AF.986.02**


Size  
**16,0cm**



**Schmieden**

Article No.  
**AF.993.19**

Size  
**19,0cm**



**Schmieden-Payr**

Article No.  
**AF.994.19**

Size  
**19,0cm**



**König 19,5cm**

Article No.  
**AF.995.03**  
**AF.995.05**  
**AF.995.07**

Size  
**2,5mm**  
**5,0mm**  
**7,0mm**



# SLITE TECH INSTRUMENTS GERMANY

## LIGATURE CONDUCTORS



**Brunner-Ferrozell**

Article No.	Size
AF.987.18	18,5cm



**Kirschner-Ferrozell**

Article No.	Size
AF.987.27	27,0cm



**Ferrozell**

Article No.	Size
AF.992.24	24,0cm

11mm breit



**MEDE Hohlsonde gebogen**

Article No.	Size
AF.990.01	19,0cm

**MEDE Hohlsonde gerade**


Article No.	Size
AF.990.02	19,0cm

AF



## KONTAKTIERE JETZT

 [info@slitetechnstruments.de](mailto:info@slitetechnstruments.de)

 +49 162 9224643

 [slitetechnstruments01@gmail.com](mailto:slitetechnstruments01@gmail.com)

 Osterkamp 10 33605,  
Bielefeld Germany

 [www.slitetechnstruments.de](http://www.slitetechnstruments.de)